

ON-SITE OXYGEN SOLUTIONS

OXYGEN STORAGE AND DISTRIBUTION

PCI offers the complete package when it comes to oxygen generation, storage, and distribution. Our oxygen solutions can help you become completely independent from the chores and costs of a cylinder re-supply chain or you can use them as a back-up system for emergencies.

STORAGE



MOST

Medical Oxygen Storage Tank

Oxygen storage and distribution system consisting of lightweight, high strength, high pressure composite wound brass lined cylinders secured within a rugged water resistant case.



SC MOST

Single Cylinder Most

Ready-to-go multi-patient 19,000 liter O₂ storage & delivery system on wheels. Provides almost 20,000 liters of oxygen at 3,000 psig service pressure.



HOBS

Hospital Oxygen Backup Supply

Designed to provide oxygen continuity in an emergency. Four-legged stand assembly with a manifold unit for grouping up to eight "H" size oxygen cylinders.

DISTRIBUTION



CODS

Civilian Oxygen Distribution System

PCI's oxygen delivery system designed for MFHs uses a modular approach so that it can be easily tailored and scaled to MFH bed capacities.



PODS

Patient Oxygen Distribution Systems

A system of flow regulators to mount at patients bedside and multiple hoses with quick disconnect fittings to provide the distribution of the oxygen at a regulated pressure and flow.



SSODS

Surgical Suite Oxygen Distribution System

Provides oxygen to the surgical suite for operating equipment such as ventilators, respirators and anesthesia machines.